



**LIFE IS TIME™**  
 When minutes count  
**EMERGENCY OXYGEN**  
 can be the difference between  
**LIFE OR DEATH**



*Only the very best and complete  
**EMERGENCY OXYGEN UNIT**  
 provides supplemental oxygen  
 for First Aid in **TWO WAYS:***

- **INHALATOR**  
 ... for a breathing victim  
*and*
- **CPR RESUSCITATOR adjunct**  
 ... for a non-breathing victim



LIFE® OxygenPac features a LIFE® mask for "mouth-to-mask" resuscitation with a one-way valve. This eliminates direct contact with the victim's mouth and nose while resuscitating a non-breathing person. Maximum hygiene is provided by the one-way valve, which eliminates transmittal of fluids and contagious diseases. LIFE® mask fits both adult and child.

(Patented)

**LIFE® OxygenPac MEETS OR EXCEEDS ALL GOVERNMENT AND FDA REQUIREMENTS: A NON-PRESCRIPTION DEVICE INTENDED FOR EMERGENCY USE.**

[www.LIFECorporation.com](http://www.LIFECorporation.com)

Patents Approved and Pending



*This unit  
 could save  
 your life.*

**LIFE® OxygenPac**  
*should be part of your first aid program.*

- SIMPLE and EASY TO USE, with external Two-Step Instructions and Clear Cover
- Wall Mounted for Recognition and Accessibility always available on display
- Easily Portable with Convenient Handle
- Simple "On-Off" Control Lever
- Constant Reading Supply Gauge  
 90 Minute Supply\*
- Resuscitation LIFE® mask with One-Way Valve Fits both adult and child
- REFILLABLE, low cost, by any compressed gas distributor with universal CGA #870 fitting
- Tight Seal Protects Unit from Dust and Moisture
- Units come Ready to Use - Full of oxygen with wall brackets and screws

MODEL	FEATURES
LIFE-6FF	6 LPM (Liters Per Minute) Fixed Flow

<b>LIFE-612</b>	6 & 12 LPM - "NORM" & "HIGH" (2 settings) provides up to recommended 100% inspired O2
-----------------	------------------------------------------------------------------------------------------

LIFE-025	0 to 25 LPM Variable Flow (for EMT's) 12 flow settings - 0, .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, 25
----------	------------------------------------------------------------------------------------------------------

Made in the U.S.A.

©2002 LIFE Corporation, Milwaukee

**IDEAL LOCATIONS FOR THE LIFE® OxygenPac  
... WHEREVER THERE ARE PEOPLE!**



**INDUSTRIAL/FACTORY**



**PUBLIC FACILITIES**



**ENTERTAINMENT/RECREATION**



**OFFICES**

LIFE® OxygenPac is attractive. Prominently displayed, it makes the simple but strong statement, "This Organization Cares For Its Employees And Visitors."

**American Medical Association  
GUIDELINES ON OXYGEN EQUIPMENT IN CPR ...**  
as published by the American Heart Association

- "... supplemental oxygen should be used as soon as it is available ..."
- "... oxygen should not be withheld for fear of suppressing respiration ..."
- "The highest possible oxygen concentration (preferable 100%) should be administered as soon as possible to all patients with cardiac or pulmonary arrest or other patients with suspected hypoxemia regardless of cause."
- "Masks equipped with one-way valve ... are recommended for mouth to mask ventilation."



**EMERGENCY OXYGEN  
INSTRUCTIONS**

1. Move lever → ON
2. Place mask over nose and mouth



IF BREATHING HAS STOPPED, START MOUTH-TO-MASK RESUSCITATION (BLOW UNTIL CHEST RISES, EVERY 5 SECONDS)

**612 MODEL LIFE® OxygenPac**

- SAME SUPERIOR FEATURES of standard LIFE® OxygenPac
- PLUS: NORM & HIGH FLOW RATES  
Regulator with "Click-Stop" two settings
  - "NORM" - 6 LITERS PER MINUTE for GENERAL emergency situations delivers approximately 45% inspiration concentration to the average adult breathing victim as an INHALATOR
  - "HIGH" - 12 LITERS PER MINUTE for EXTREME emergency situations delivers oxygen concentrations approaching the recommended 100% inspiration concentration as an INHALATOR for breathing victims and as a RESUSCITATOR for oxygen supplemented mouth-to-mask resuscitation for a non-breathing victim.
- SIMPLE OPERATION IS EASY-TO-USE  
LIFE® OxygenPac 612 MODEL cannot be turned to below 6 LPM (Lower flow rates are not advised when used with mask delivery.)

**5 YEAR WARRANTY**

Our warranty includes provisions that all parts and materials are warranted against failure due to defective material or poor workmanship for FIVE YEARS from date of purchase.

**LIFE® OxygenPac TECHNICAL SPECIFICATIONS**

<p><b>UL</b> APPROVED REGULATOR</p> <p>Flow length.....90 minutes ± U.S.P oxygen Liter flow.....6 LPM constant Style.....Piston</p>		<p><b>CYLINDER (included)</b> <b>566 liters</b></p> <p>Fill capacity.....2015 PSI "Fill" pressure.....3358 PSI Test pressure.....3358 PSI Material.....REFILLABLE seamless steel</p>
<p><b>OXYGEN PAC UNIT</b></p> <p>H x W x D.....<b>20" x 6 3/8" x 9"</b> Weight.....<b>16 lbs.</b> Component material: Case.....Polyethylene Viewing cover.....Lexan® polycarbonate Tight seal protects unit from dust and moisture</p>		<p><b>Dimensions:</b> Outside diameter.....5 1/4" Overall length.....14" Minimum wall thickness.....0.108"</p>
<p><b>GAUGE</b></p> <p>Readout.....Constant Cylinder contents "Full" indicator.....Green wedge Maximum pressure rating...4000 PSI</p>		<p><b>MASK</b></p> <p>Style.....Mouth-to-mask resuscitator with one-way flow valve Material.....Clear vinyl</p>
<p><b>SHUT-OFF VALVE</b></p> <p>"On-Off" lever.....Adjustable ratchet Thread.....1/2" - 14 NGT taper Safety burst relief value.....3025 PSI Refill adaptor.....C.G.A. #870</p>		<p><b>TUBING</b></p> <p>Length.....7 ft. Crimp-proof star diameter</p> <p>Brackets and Fasteners for Wall Mounting Included with Every Unit</p> <p>REFILLABLE by Any Medical Industrial Gas Distributor with Universally Standard CGA #870 Fitting</p>

**ORDERING INFORMATION  
MODEL # FEATURES**

LIFE-6FF	6 LPM (Liters Per Minute) Fixed Flow
<b>LIFE-612</b>	<b>6 &amp; 12 LPM - "NORM" &amp; "HIGH" (2 settings)</b> provides up to recommended 100% inspired O2
LIFE-025	0 to 25 LPM Variable Flow (for EMT's) 12 flow settings - 0, .25, .5, 1, 2, 3, 4, 6, 8, 10, 15, 25
LIFE-101	Spare Cylinder Assembly (cylinder, valve, gauge, O2)
LIFE-102	LIFE® mask & 1-way valve Fits both adult and child (mask, valve, strap, O2 npl)
LIFE-BAG	LIFE® Bag/mask manual Resuscitator with O2 inlet

(All units complete with Cylinder, Oxygen, Brackets and Fasteners for Wall Mounting)

Uncompromisingly manufactured with quality materials, technology and workmanship



Incontestable trademarks of LIFE Corporation



PATENTS APPROVED AND PENDING Made in the U.S.A.  
© 2002 LIFE Corporation, Milwaukee

**TO ORDER CALL:**

"World's Largest AED Source"

**AEDSuperstore®**  
an Allied 100™ company

**800.544.0048**  
**help@aeds.com**



**www.AEDSuperstore.com**

1800 US Hwy 51 N • Woodruff, WI 54568 • 715.358.2329